



AVPR1A Antibody, HRP conjugated

Product Code	CSB-PA002468LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37288
Immunogen	Recombinant Human Vasopressin V1a receptor protein (1-52AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vasopressin V1a receptor (V1aR) (AVPR V1a) (Antidiuretic hormone receptor 1a) (Vascular/hepatic-type arginine vasopressin receptor), AVPR1A, AVPR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	AVPR1A