



CDS2 Antibody, HRP conjugated

Product Code	CSB-PA005123LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95674
Immunogen	Recombinant Human Phosphatidate cytidyltransferase 2 protein (354-445AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidate cytidyltransferase 2 (EC 2.7.7.41) (CDP-DAG synthase 2) (CDP-DG synthase 2) (CDP-diacylglycerol synthase 2) (CDS 2) (CDP-diglyceride pyrophosphorylase 2) (CDP-diglyceride synthase 2) (CTP:phosphatidate cytidyltransferase 2), CDS2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CDS2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.