



ETV5 Antibody, Biotin conjugated

Product Code	CSB-PA007859LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41161
Immunogen	Recombinant Human ETS translocation variant 5 protein (124-265AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ETS translocation variant 5 (Ets-related protein ERM), ETV5, ERM
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ETV5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.