



F13B Antibody, Biotin conjugated

Product Code	CSB-PA007922LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05160
Immunogen	Recombinant Human Coagulation factor XIII B chain protein (260-403AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coagulation factor XIII B chain (Fibrin-stabilizing factor B subunit) (Protein-glutamine gamma-glutamyltransferase B chain) (Transglutaminase B chain), F13B
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	F13B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.