



FBP2 Antibody, FITC conjugated

Product Code	CSB-PA008460LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00757
Immunogen	Recombinant Human Fructose-1,6-bisphosphatase isozyme 2 protein (283-339AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fructose-1,6-bisphosphatase isozyme 2 (FBPase 2) (EC 3.1.3.11) (D-fructose-1,6-bisphosphate 1-phosphohydrolase 2) (Muscle FBPase), FBP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FBP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.