



FCGR3B Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA008544LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O75015 |
| Immunogen | Recombinant Human Low affinity immunoglobulin gamma Fc region receptor III-B protein (18-125AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Low affinity immunoglobulin gamma Fc region receptor III-B (Fc-gamma RIII-beta) (Fc-gamma RIII) (Fc-gamma RIIIb) (FcRIII) (FcRIIIb) (FcR-10) (IgG Fc receptor III-1) (CD antigen CD16b), FCGR3B, CD16B FCG3 FCGR3 IGFR3 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | FCGR3B |