



FCGR3B Antibody, HRP conjugated

Product Code	CSB-PA008544LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75015
Immunogen	Recombinant Human Low affinity immunoglobulin gamma Fc region receptor III-B protein (18-125AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Low affinity immunoglobulin gamma Fc region receptor III-B (Fc-gamma RIII-beta) (Fc-gamma RIII) (Fc-gamma RIIIb) (FcRIII) (FcRIIIb) (FcR-10) (IgG Fc receptor III-1) (CD antigen CD16b), FCGR3B, CD16B FCG3 FCGR3 IGFR3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FCGR3B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.