



FOLR1 Antibody, Biotin conjugated

Product Code	CSB-PA008784LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15328
Immunogen	Recombinant Human Folate receptor alpha protein (25-234AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Folate receptor alpha (FR-alpha) (Adult folate-binding protein) (FBP) (Folate receptor 1) (Folate receptor, adult) (KB cells FBP) (Ovarian tumor-associated antigen MOv18), FOLR1, FOLR
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	FOLR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.