



FOLR2 Antibody, FITC conjugated

Product Code	CSB-PA008786LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14207
Immunogen	Recombinant Human Folate receptor beta protein (35-167AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Folate receptor beta (FR-beta) (Folate receptor 2) (Folate receptor, fetal/placental) (Placental folate-binding protein) (FBP), FOLR2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FOLR2