



GNA12 Antibody, FITC conjugated

Product Code	CSB-PA009584LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03113
Immunogen	Recombinant Human Guanine nucleotide-binding protein subunit alpha-12 protein (112-270AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanine nucleotide-binding protein subunit alpha-12 (G alpha-12) (G-protein subunit alpha-12), GNA12
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GNA12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.