



GPHN Antibody, HRP conjugated

Product Code	CSB-PA009716LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQX3
Immunogen	Recombinant Human Gephyrin protein (360-460AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gephyrin [Includes: Molybdopterin adenylyltransferase (MPT adenylyltransferase) (EC 2.7.7.75) (Domain G); Molybdopterin molybdenumtransferase (MPT Mo-transferase) (EC 2.10.1.1) (Domain E)], GPHN, GPH KIAA1385
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GPHN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.