



IGKC Antibody, Biotin conjugated

Product Code	CSB-PA011340LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01834
Immunogen	Recombinant Human Immunoglobulin kappa constant protein (1-107AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Immunoglobulin kappa constant (Ig kappa chain C region) (Ig kappa chain C region AG) (Ig kappa chain C region CUM) (Ig kappa chain C region EU) (Ig kappa chain C region OU) (Ig kappa chain C region ROY) (Ig kappa chain C region TI), IGKC
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IGKC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.