



ELP1 Antibody, Biotin conjugated

Product Code	CSB-PA011571LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95163
Immunogen	Recombinant Human Elongator complex protein 1 protein (1107-1258AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Elongator complex protein 1 (ELP1) (IkappaB kinase complex-associated protein) (IKK complex-associated protein) (p150), ELP1, IKAP, IKBKAP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ELP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.