



MAP4K2 Antibody, Biotin conjugated

Product Code	CSB-PA013437LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12851
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase kinase 2 protein (386-469AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase kinase 2 (EC 2.7.11.1) (B lymphocyte serine/threonine-protein kinase) (Germinal center kinase) (GC kinase) (MAPK/ERK kinase kinase kinase 2) (MEK kinase kinase 2) (MEKKK 2) (Rab8-interacting protein), MAP4K2, GCK RAB8IP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP4K2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.