



MCC Antibody, FITC conjugated

Product Code	CSB-PA013571LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23508
Immunogen	Recombinant Human Colorectal mutant cancer protein (479-611AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Colorectal mutant cancer protein (Protein MCC), MCC
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	MCC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.