



NLRP6 Antibody, HRP conjugated

Product Code	CSB-PA015874LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P59044
Immunogen	Recombinant Human NACHT, LRR and PYD domains-containing protein 6 protein (154-306AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NACHT, LRR and PYD domains-containing protein 6 (Angiotensin II/vasopressin receptor) (PYRIN-containing APAF1-like protein 5), NLRP6, NALP6 PYPAF5
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NLRP6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.