



PEX11A Antibody, FITC conjugated

Product Code	CSB-PA017795LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75192
Immunogen	Recombinant Human Peroxisomal membrane protein 11A protein (106-219AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisomal membrane protein 11A (HsPEX11p) (28 kDa peroxisomal integral membrane protein) (PMP28) (Peroxin-11A) (Peroxisomal biogenesis factor 11A) (Protein PEX11 homolog alpha) (PEX11-alpha), PEX11A, PEX11
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PEX11A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.