




PGLYRP3 Antibody

Product Code	CSB-PA017864LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LB9
Immunogen	Recombinant Human Peptidoglycan recognition protein 3 protein (90-222AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidoglycan recognition protein 3 (Peptidoglycan recognition protein I-alpha) (PGLYRPIalpha) (PGRP-I-alpha) (Peptidoglycan recognition protein intermediate alpha), PGLYRP3, PGRPIA
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PGLYRP3

Image



Western Blot

Positive WB detected in: MCF-7 whole cell lysate, Colo320 whole cell lysate

All lanes: PGLYRP3 antibody at 4.86µg/ml
Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 38 kDa

Observed band size: 38 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.