



PGLYRP3 Antibody, FITC conjugated

Product Code	CSB-PA017864LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LB9
Immunogen	Recombinant Human Peptidoglycan recognition protein 3 protein (90-222AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidoglycan recognition protein 3 (Peptidoglycan recognition protein I-alpha) (PGLYRP1alpha) (PGRP-I-alpha) (Peptidoglycan recognition protein intermediate alpha), PGLYRP3, PGRPIA
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PGLYRP3