



PGLYRP3 Antibody, Biotin conjugated

Product Code	CSB-PA017864LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LB9
Immunogen	Recombinant Human Peptidoglycan recognition protein 3 protein (90-222AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidoglycan recognition protein 3 (Peptidoglycan recognition protein I-alpha) (PGLYRP1alpha) (PGRP-I-alpha) (Peptidoglycan recognition protein intermediate alpha), PGLYRP3, PGRP1A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PGLYRP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.