



RAB7B Antibody, FITC conjugated

Product Code	CSB-PA019220LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AH8
Immunogen	Recombinant Human Ras-related protein Rab-7b protein (100-199AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-7b, RAB7B
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	RAB7B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.