



REM1 Antibody, Biotin conjugated

Product Code	CSB-PA019559LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75628
Immunogen	Recombinant Human GTP-binding protein REM 1 protein (2-153AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GTP-binding protein REM 1 (GTPase-regulating endothelial cell sprouting) (Rad and Gem-like GTP-binding protein 1), REM1, GES REM
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	REM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.