



RIN3 Antibody, FITC conjugated

Product Code	CSB-PA019728LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TB24
Immunogen	Recombinant Human Ras and Rab interactor 3 protein (245-465AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras and Rab interactor 3 (Ras interaction/interference protein 3), RIN3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RIN3