



RPS23 Antibody, Biotin conjugated

Product Code	CSB-PA020400LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62266
Immunogen	Recombinant Human 40S ribosomal protein S23 protein (2-143AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S23 (Small ribosomal subunit protein uS12), RPS23
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPS23
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.