



SASH1 Antibody, HRP conjugated

Product Code	CSB-PA020714LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94885
Immunogen	Recombinant Human SAM and SH3 domain-containing protein 1 protein (76-191AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SAM and SH3 domain-containing protein 1 (Proline-glutamate repeat-containing protein), SASH1, KIAA0790 PEPE1
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	SASH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.