



SEC14L2 Antibody, HRP conjugated

Product Code	CSB-PA020937LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76054
Immunogen	Recombinant Human SEC14-like protein 2 protein (233-358AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SEC14-like protein 2 (Alpha-tocopherol-associated protein) (TAP) (hTAP) (Squalene transfer protein) (Supernatant protein factor) (SPF), SEC14L2, C22orf6 KIAA1186 KIAA1658
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SEC14L2