



SEC14L2 Antibody, Biotin conjugated

Product Code	CSB-PA020937LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76054
Immunogen	Recombinant Human SEC14-like protein 2 protein (233-358AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SEC14-like protein 2 (Alpha-tocopherol-associated protein) (TAP) (hTAP) (Squalene transfer protein) (Supernatant protein factor) (SPF), SEC14L2, C22orf6 KIAA1186 KIAA1658
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SEC14L2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.