



SLC12A3 Antibody, Biotin conjugated

Product Code	CSB-PA021384LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55017
Immunogen	Recombinant Human Solute carrier family 12 member 3 protein (791-952AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 12 member 3 (Na-Cl cotransporter) (NCC) (Na-Cl symporter) (Thiazide-sensitive sodium-chloride cotransporter), SLC12A3, NCC TSC
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC12A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.