



SLC22A5 Antibody, HRP conjugated

Product Code	CSB-PA021460LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76082
Immunogen	Recombinant Human Solute carrier family 22 member 5 protein (42-142AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 22 member 5 (High-affinity sodium-dependent carnitine cotransporter) (Organic cation/carnitine transporter 2), SLC22A5, OCTN2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC22A5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.