




SLCO1A2 Antibody

Product Code	CSB-PA021741LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46721
Immunogen	Recombinant Human Solute carrier organic anion transporter family member 1A2 protein (406-513AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier organic anion transporter family member 1A2 (OATP-A) (Organic anion-transporting polypeptide 1) (OATP-1) (Sodium-independent organic anion transporter) (Solute carrier family 21 member 3), SLCO1A2, OATP OATP1 OATP1A2 SLC21A3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLCO1A2

Image



Immunofluorescence staining of SH-SY5Y cells with CSB-PA021741LA01HU at 1:166, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.