



SLCO1A2 Antibody, HRP conjugated

Product Code	CSB-PA021741LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46721
Immunogen	Recombinant Human Solute carrier organic anion transporter family member 1A2 protein (406-513AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier organic anion transporter family member 1A2 (OATP-A) (Organic anion-transporting polypeptide 1) (OATP-1) (Sodium-independent organic anion transporter) (Solute carrier family 21 member 3), SLCO1A2, OATP OATP1 OATP1A2 SLC21A3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLCO1A2