



SNRPE Antibody, HRP conjugated

Product Code	CSB-PA022336LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62304
Immunogen	Recombinant Human Small nuclear ribonucleoprotein E protein (1-92AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Small nuclear ribonucleoprotein E (snRNP-E) (Sm protein E) (Sm-E) (SmE), SNRPE
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SNRPE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.