



SPAG6 Antibody, HRP conjugated

Product Code	CSB-PA022467LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75602
Immunogen	Recombinant Human Sperm-associated antigen 6 protein (343-428AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sperm-associated antigen 6 (Protein PF16 homolog) (Repro-SA-1) (Sperm flagellar protein), SPAG6, PF16
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SPAG6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.