



SSR2 Antibody, FITC conjugated

Product Code	CSB-PA022718LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43308
Immunogen	Recombinant Human Translocon-associated protein subunit beta protein (18-149AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Translocon-associated protein subunit beta (TRAP-beta) (Signal sequence receptor subunit beta) (SSR-beta), SSR2, TRAPB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SSR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.