



SUCLG1 Antibody, HRP conjugated

Product Code	CSB-PA022919LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53597
Immunogen	Recombinant Human Succinate--CoA ligase [ADP/GDP-forming] subunit alpha, mitochondrial protein (210-316AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Succinate--CoA ligase [ADP/GDP-forming] subunit alpha, mitochondrial (EC 6.2.1.4) (EC 6.2.1.5) (Succinyl-CoA synthetase subunit alpha) (SCS-alpha), SUCLG1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SUCLG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.