



SV2C Antibody, FITC conjugated

Product Code	CSB-PA022980LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q496J9
Immunogen	Recombinant Human Synaptic vesicle glycoprotein 2C protein (459-524AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Synaptic vesicle glycoprotein 2C, SV2C, KIAA1054
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SV2C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.