



DCSTAMP Antibody, Biotin conjugated

Product Code	CSB-PA023625LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H295
Immunogen	Recombinant Human Dendritic cell-specific transmembrane protein (119-209AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dendritic cell-specific transmembrane protein (DC-STAMP) (hDC-STAMP) (Dendrocyte-expressed seven transmembrane protein) (IL-four-induced protein) (FIND) (Transmembrane 7 superfamily member 4), DCSTAMP, TM7SF4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DCSTAMP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.