



TMEM132D Antibody, Biotin conjugated

Product Code	CSB-PA023702LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14C87
Immunogen	Recombinant Human Transmembrane protein 132D protein (796-909AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 132D (Mature oligodendrocytes transmembrane protein) (Mature OL transmembrane protein), TMEM132D, HBE120 KIAA1944 MOLT
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TMEM132D
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.