



TNFRSF18 Antibody, Biotin conjugated

Product Code	CSB-PA023975LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5U5
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 18 protein (90-162AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 18 (Activation-inducible TNFR family receptor) (Glucocorticoid-induced TNFR-related protein) (CD antigen CD357), TNFRSF18, AITR GITR
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TNFRSF18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.