



TNFRSF18 Antibody, Biotin conjugated

Product Code	CSB-PA023975LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5U5
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 18 protein (90-162AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 18 (Activation-inducible TNFR family receptor) (Glucocorticoid-induced TNFR-related protein) (CD antigen CD357), TNFRSF18, AITR GITR
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TNFRSF18