



TTC26 Antibody, Biotin conjugated

Product Code	CSB-PA025219LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0AVF1
Immunogen	Recombinant Human Intraflagellar transport protein 56 protein (5-122AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intraflagellar transport protein 56 (Tetratricopeptide repeat protein 26) (TPR repeat protein 26), TTC26, IFT56
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TTC26
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.