



ZNF3 Antibody, HRP conjugated

Product Code	CSB-PA026651LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17036
Immunogen	Recombinant Human Zinc finger protein 3 protein (92-210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 3 (Zinc finger protein HF.12) (Zinc finger protein HZF3.1) (Zinc finger protein KOX25), ZNF3, KOX25
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.