



VTN Antibody, Biotin conjugated

Product Code	CSB-PA029544ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04004
Immunogen	Recombinant Human Vitronectin protein (1-115AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vitronectin (VN) (S-protein) (Serum-spreading factor) (V75) [Cleaved into: Vitronectin V65 subunit; Vitronectin V10 subunit; Somatomedin-B], VTN
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	VTN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.