



TFR2 Antibody, FITC conjugated

Product Code	CSB-PA14729C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UP52
Immunogen	Recombinant Human Transferrin receptor protein 2 protein (534-605AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transferrin receptor protein 2 (TfR2), TFR2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TFR2