



FAN1 Antibody, HRP conjugated

Product Code	CSB-PA22859B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2M0
Immunogen	Recombinant Human Fanconi-associated nuclease 1 protein (167-373AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fanconi-associated nuclease 1 (EC 3.1.21.-) (EC 3.1.4.1) (FANCD2/FANCI-associated nuclease 1) (hFAN1) (Myotubularin-related protein 15), FAN1, KIAA1018 MTMR15
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FAN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.