



FAN1 Antibody, HRP conjugated

Product Code	CSB-PA22859B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2M0
Immunogen	Recombinant Human Fanconi-associated nuclease 1 protein (167-373AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fanconi-associated nuclease 1 (EC 3.1.21.-) (EC 3.1.4.1) (FANCD2/FANCI-associated nuclease 1) (hFAN1) (Myotubularin-related protein 15), FAN1, KIAA1018 MTMR15
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FAN1