



ITIH3 Antibody, HRP conjugated

Product Code	CSB-PA22869B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06033
Immunogen	Recombinant Human Inter-alpha-trypsin inhibitor heavy chain H3 protein (512-651AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inter-alpha-trypsin inhibitor heavy chain H3 (ITI heavy chain H3) (ITI-HC3) (Inter-alpha-inhibitor heavy chain 3) (Serum-derived hyaluronan-associated protein) (SHAP), ITIH3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ITIH3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.