



KLHL20 Antibody, Biotin conjugated

Product Code	CSB-PA22879D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2M5
Immunogen	Recombinant Human Kelch-like protein 20 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kelch-like protein 20 (Kelch-like ECT2-interacting protein) (Kelch-like protein X), KLHL20, KLEIP KLHLX
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	KLHL20
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.