



VKORC1 Antibody, HRP conjugated

Product Code	CSB-PA22899B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BQB6
Immunogen	Recombinant Human Vitamin K epoxide reductase complex subunit 1 protein (30-74AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vitamin K epoxide reductase complex subunit 1 (EC 1.17.4.4) (Vitamin K1 2,3-epoxide reductase subunit 1), VKORC1, VKOR
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	VKORC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.