



VKORC1 Antibody, Biotin conjugated

Product Code	CSB-PA22899D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BQB6
Immunogen	Recombinant Human Vitamin K epoxide reductase complex subunit 1 protein (30-74AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vitamin K epoxide reductase complex subunit 1 (EC 1.17.4.4) (Vitamin K1 2,3-epoxide reductase subunit 1), VKORC1, VKOR
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	VKORC1